


AL300RV


RECREATIONAL VEHICLE LIFT



Easy, one switch operation

4 self-tensioning Q-Strait retractors

Rear wheel cradle



REVOLUTIONARY BY DESIGN

Engineered to handle the stress of wind and heavy travel with robust hitch adapter to manage higher forces. Larger tubing in the platform avoids structural fatigue in off-road conditions.

AL300RV PRODUCT SPECIFICATIONS

LIFTING CAPACITY 350 lb

ACCOMMODATES - Power wheelchairs
- Scooters

HITCH CLASS III or greater

POWERED - Manual fold
- Power lifting

AVAILABLE OPTION(S) - Battery pack
- Cover

INSTALLED WEIGHT 82 lb

PLATFORM SIZE 26" x 39" adjustable cradle extends up to 46"

WARRANTY 3-year

INCLUDED - AL107 hitch post
- H300102 class III reinforced hitch adapter
- License plate mount

SECUREMENT - 4 self-tensioning Q'Strait retractors
- Rear wheel cradle



Harmar.com
800.833.0478

25SEP2020 | MKT-000129 REV B

