

AL001

UNIVERSAL MANUAL CARRIER



AFFORDABLE, FULLY-MANUAL SOLUTION

Economical, fully manual solution carries virtually any power wheelchair or scooter.

AL001 PRODUCT SPECIFICATIONS

LIFTING CAPACITY 500 lb

PLATFORM SIZE 28" x 48" platform with a 60" folding ramp

HITCH CLASS II or III

AVAILABLE OPTION(S) - Cover
- Tie-down straps

INSTALLED WEIGHT 84 lb

ACCOMMODATES - Manual wheelchairs
- Power wheelchairs
- Scooters

WARRANTY 1-year



Harmar.com
800.833.0478

08AUG2020 | MKT-000129 REV B

