

## AL015

### UNIVERSAL MICRO SCOOTER LIFT



#### **ADJUSTABLE, LIGHTWEIGHT, PROFILE PLATFORM THAT MAXIMIZES VISIBILITY**

Ideal for compact three or four wheel scooters and is compatible with many smaller vehicles.

Adjustable cradles are able to accommodate a wide variety of scooters up to 22" wide and wheelbases from 26" to 36".



Adjustable  
hold-down arm

Optional  
4-wheel cradle  
(AL0154)



The AL105 is a great solution for my little scooter and car. I have great visibility behind when its folded up and it's light enough that I don't need a truck.

## AL015 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	150 lb	<b>ACCOMMODATES</b>	Smaller scooters up to 36" wheelbase
<b>INSTALLED WEIGHT</b>	- 65 lb (3-wheel cradle standard) - 70 lb (4-wheel cradle option)	<b>PLATFORM SIZE</b>	22.5" x 40"
<b>HITCH CLASS</b>	I or greater	<b>POWERED</b>	Fully powered
<b>WARRANTY</b>	3-year	<b>INCLUDED</b>	AL126 class II adaptor
<b>AVAILABLE OPTION(S)</b>	- Battery pack - Swing-away - Scooter cover - 4-wheel cradle	<b>SECUREMENT</b>	- Automatic hold-down foot - Front and rear wheel cradle