


AL560 / AL560XL

AUTOMATIC UNIVERSAL POWER WHEELCHAIR LIFT



One switch operation and key lock

Adjustable wheel cradles and platform skins

Hold-down arm

ADJUSTABLE POWER CHAIR SOLUTION

The easy way to secure a wide range of popular power chairs. Hold-Down arm is padded to hold securely while cushioning chair. The 26.25" x 38.5" platform (AL560) and 30.75" x 38.5" (AL560XL) platform adjusts for a custom fit. Manual crank backup and license plate mount are included.

Swing-away option available.

AL560 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	Power wheelchairs
INSTALLED WEIGHT	89 lb
PLATFORM SIZE	26.25" x 38.5" adjustable platform
HITCH CLASS	II or greater
POWERED	Fully powered
WARRANTY	3-year
INCLUDED	<ul style="list-style-type: none">- AL107 hitch post- AL123 class III hitch adapter- License plate mount
SECUREMENT	<ul style="list-style-type: none">- Padded hold-down arm- Wheel cradles
AVAILABLE OPTION(S)	<ul style="list-style-type: none">- Battery pack- Cover- Swing-away

AL560XL PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	Larger power wheelchairs
INSTALLED WEIGHT	101 lb
PLATFORM SIZE	30.75" x 38.5" adjustable platform
HITCH CLASS	III or greater
POWERED	Fully powered
WARRANTY	3-year
INCLUDED	<ul style="list-style-type: none">- AL107 hitch post- AL123HD class III heavy duty hitch adapter- License plate mount
SECUREMENT	<ul style="list-style-type: none">- Padded hold-down arm- Wheel cradles
AVAILABLE OPTION(S)	<ul style="list-style-type: none">- Battery pack- Cover- Swing-away



Harmar.com
800.833.0478

17SEP2020 | MKT-000176 REV A

