

# AL580 / AL580XL

## NEXT GENERATION POWER WHEELCHAIR LIFT



One switch operation and key lock

Hold-down arm

Wheel cradle

### LIGHTEST SOLUTION ON THE MARKET

The AL580 and AL580XL are compact, durable outside vehicle platform lifts that safely and securely transports power chairs.

The AL580 series is easy to operate and built to last for years of trouble-free use. Designed to meet the needs of consumers who need a quick, automatic securing option without making modifications to their power wheelchair, the stylish and compact detailing of this lift makes it a must have for power wheelchair owners.

### AL580 PRODUCT SPECIFICATIONS

|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>LIFTING CAPACITY</b>    | 350 lb                                                                                                                                   |
| <b>ACCOMMODATES</b>        | Mid-wheel drive power wheelchairs with 10" drive wheels                                                                                  |
| <b>INSTALLED WEIGHT</b>    | 63 lb                                                                                                                                    |
| <b>PLATFORM SIZE</b>       | 16.5" x 27"                                                                                                                              |
| <b>HITCH CLASS</b>         | II or greater                                                                                                                            |
| <b>POWERED</b>             | Fully powered                                                                                                                            |
| <b>WARRANTY</b>            | 3-year                                                                                                                                   |
| <b>INCLUDED</b>            | <ul style="list-style-type: none"><li>- AL107 hitch post</li><li>- AL123 class III hitch adapter</li><li>- License plate mount</li></ul> |
| <b>SECUREMENT</b>          | <ul style="list-style-type: none"><li>- Padded hold-down arm</li><li>- Wheel cradle</li></ul>                                            |
| <b>AVAILABLE OPTION(S)</b> | <ul style="list-style-type: none"><li>- Battery pack</li><li>- Cover</li><li>- License plate riser</li><li>- Swing-away</li></ul>        |

### AL580XL PRODUCT SPECIFICATIONS

|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>LIFTING CAPACITY</b>    | 350 lb                                                                                                                                   |
| <b>ACCOMMODATES</b>        | Mid-wheel drive power wheelchairs with 14" drive wheels                                                                                  |
| <b>INSTALLED WEIGHT</b>    | 63 lb                                                                                                                                    |
| <b>PLATFORM SIZE</b>       | 17.5" x 27"                                                                                                                              |
| <b>HITCH CLASS</b>         | II or greater                                                                                                                            |
| <b>POWERED</b>             | Fully powered                                                                                                                            |
| <b>WARRANTY</b>            | 3-year                                                                                                                                   |
| <b>INCLUDED</b>            | <ul style="list-style-type: none"><li>- AL107 hitch post</li><li>- AL123 class III hitch adapter</li><li>- License plate mount</li></ul> |
| <b>SECUREMENT</b>          | <ul style="list-style-type: none"><li>- Padded hold-down arm</li><li>- Wheel cradle</li></ul>                                            |
| <b>AVAILABLE OPTION(S)</b> | <ul style="list-style-type: none"><li>- Battery pack</li><li>- Cover</li><li>- License plate riser</li><li>- Swing-away</li></ul>        |



Harmar.com  
800.833.0478

17SEP2020 | MKT-000177 REV A

