## **PINNACLE SSL-SL600**

## EVALUATION/ORDER FORM

Quote Order	
Date:	Account Number:
Sales Representative:	Projected Close Date:
Dealer Name:	Purchase Order Number:
Address:	Ship To Address:
City: State: Zip:	
Phone: Fax:	Contact:
Email:	Phone:
Shipping Method: Pick up Common Carrier	Opportunity Reference / Project Name:
PRIMARY MODEL	
QTY: SSL-SL600: Standard 15'6" rail	
CONFIGURATION	
QTY: SL600-PSPF: Power Swivel + Power Footrest bundle	QTY: SSL-SL600FR-LH: Left hand folding rail; 15' 6" rail
QTY: SL600-PS: Power Swivel	QTY: SSL-SL600FR-RH: Right hand folding rail; 15' 6" rail
QTY: SL600-PF: Power Footrest	
TRACK, PARTS AND OPTIONS	
TRACK	PARTS AND OPTIONS
QTY: SSC-KNT: No track or brackets in lieu of std	QTY: SSL-KEY: Key lock on lift (chassis)
QTY: SSL-8 TRACK: 7'9" track in lieu of std	QTY: SSL-SB1: Additional short brackets
QTY: SSL-T3-SV: 3' mid-section assembly with brackets	QTY: SSL-TB1: Additional tall brackets
QTY: SSL-T5-SV: 5' mid-section assembly with brackets	QTY: SSL-CORD-8: Charger cord extension
QTY: SSL-T8-SV: 7'9" mid-section assembly with brackets	QTY: SL600-MLT: Manual lowering tool
SSL-C: Custom cut*feet in length (including overrun)	QTY: SSL-WA: Lock out switch (WA)
* Single unit orders only Must include additional mid-section(s) above if total above 15'6".	QTY: SSL-EXT: Extension for retractable seatbelt
HOW TO MEASURE	
	MEASUREMENTS
	A to B: Inches from top landing to bottom floor
RISE	A to H: Inches from top landing to nearest obstruction
TREAD	(i.e. a door)  Rise in. Tread in. Stair width in.
	7" overrun recommended if top landing is clear
	(This is added to the A to B)
A to B	<u>Note:</u> You must order a Tall Bracket if <b>A to H is &lt; 13</b> ".
Total Rail Length	

